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
版本号：PNMK-4-KYTD012-A/0

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1. Accurately weigh 100 mg ($\pm 1\%$) of sample in 2 mL EP tube, and add 0.6 mL methanol, vortex for 30 seconds;
2. Add 100 mg glass beads and put the samples into TissueLysis II tissue grinding machine. Grind them at 25 Hz for 60 s;
3. Ultrasound at room temperature for 15 minutes;
4. Centrifuge at 25 °C for 10 min at 12 000 rpm, and the supernatant was filtered through 0.22 μ m membrane to obtain the prepared samples for LC-MS;
5. Take 20 μ L from each sample to the quality control (QC) samples (Fig. 1)*;
(These QC samples were used to monitor deviations of the analytical results from these pool mixtures and compare them to the errors caused by the analytical instrument itself)
6. Use the rest of the samples for LC-MS detection.

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Chromatographic separation was accomplished in an Thermo Vanquish system equipped with an ACQUITY UPLC® BEH C18 (100×2.1 mm, 1.7 μm, Waters) column maintained at 40 °C. The temperature of the autosampler was 8 °C. Gradient elution of analytes was carried out with 0.1% formic acid in water (A2) and 0.1% formic acid in acetonitrile (B2) at a flow rate of 0.25 mL/min. Injection of 2 μL of each sample was done after equilibration. An increasing linear gradient of solvent B (v/v) was used as follows: 0~1 min, 20% B2; 1~9 min, 20%~50% B2; 9~12 min, 50%~98% B2; 12~13.5 min, 98% B2; 13.5~14 min, 98%~20% B2; 14~20 min, 20% B2-positive model (14~17 min, 20% B2-negative model). The scanning range of UV detection wavelength is 200~400nm.

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The ESI-MSn experiments were executed on the Thermo Q Exactive HF-X mass spectrometer with the spray voltage of 3.8 kV and -2.5 kV in positive and negative modes, respectively. Sheath gas and auxiliary gas were set at 30 and 10 arbitrary units, respectively. The capillary temperature was 325 °C. The analyzer scanned over a mass range of m/z 81-1 000 for full scan at a mass resolution of 60 000. Data dependent acquisition (DDA) MS/MS experiments were performed with HCD scan. The normalized collision energy was 30 eV, 50 eV and 60 eV. Dynamic exclusion was implemented to remove some unnecessary information in MS/MS spectra.

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1. De Vos RC1, Moco S, et al. Untargeted large-scale plant metabolomics using liquid chromatography coupled to mass spectrometry. Nature Protocols (4):778-91 (2007).
2. Sangster, Timothy. et al. A pragmatic and readily implemented quality control strategy for HPLC-MS and GC-MS-based metabonomic analysis. Analyst 131.10 (2006): 1075-1078.
3. The HUSERMET Project. Analysis of Serum Samples by UPLC-MS. 2008.
4. Want, E. J. et al. Global metabolic profiling procedures for urine using UPLC-MS. Nature Protocols 5.6 (2010).

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