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1. Transfer 5, 10, 50, 100, 200, 300 μL of NIST serum into centrifuge tubes. Add ddH₂O, mixed internal standard solution and methanol according to Table. 4;
2. Vortex the six NIST serum, respectively. Centrifuge at 4 $^{\circ}\text{C}$ for 10 min at 12 000 rpm, and then transfer all the supernatant from each sample into another 2 mL centrifuge tube. Samples were concentrated to dry in vacuum. Dissolve samples with 150 μL of 80% methanol solution (-20 $^{\circ}\text{C}$), and centrifuge at 4 $^{\circ}\text{C}$ for 10 min at 12 000 rpm again to obtain the supernatant for LC-MS.

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1. Thaw all samples at 4 °C 【The insufficient samples are reduced to an equal scale】；
2. Transfer 100 μ L of each sample into 2 mL centrifuge tubes;
3. Add 100 μ L of mixed internal standard solution and 400 μ L of methanol (pre-cooled at -20 °C), then vortex for 60 s;
4. Centrifuge at 4 °C for 10 min at 12000 rpm, and then transfer 500 μ L of the supernatant from each sample into another 2 mL centrifuge tube. Samples were concentrated to dry in vacuum;
5. Dissolve samples with 150 μ L of 80% methanol solution, and centrifuge at 4 °C for 10 min at 12000 rpm again to obtain the supernatant for LC-MS;
6. Take 20 μ L from each sample to the quality control (QC) samples (Fig. 1)*;
(These QC samples were used to monitor deviations of the analytical results from these pool mixtures and compare them to the errors caused by the analytical instrument itself)
7. Use the rest of the samples for LC-MS detection.

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Chromatographic separation was accomplished in an Thermo Vanquish system equipped with an ACQUITY UPLC® HSS T3 (150 \times 2.1 mm, 1.8 μm , Waters) column maintained at 40 $^{\circ}\text{C}$. The temperature of the autosampler was 8 $^{\circ}\text{C}$. Gradient elution of analytes was carried out with 0.1% formic acid in water (A1) and 0.1% formic acid in acetonitrile (B1) or 5 mM ammonium formate in water (A3) and acetonitrile (B3) at a flow rate of 0.25 mL/min. Injection of 2 μL of each sample was done after equilibration. An increasing linear gradient of solvent A (v/v) was used as follows: 0~1 min, 2% B1/B3; 1~9 min, 2%~50% B1/B3; 9~12 min, 50%~98% B1/B3; 12~13.5 min, 98% B1/B3; 13.5~14 min, 98%~2% B1/B3; 14~20 min, 2% B1-positive model (14~17 min, 2% B3-negative model).

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The ESI-MS_n experiments were executed on the Thermo Q Exactive mass spectrometer with the spray voltage of 3.8 kV and -2.5 kV in positive and negative modes, respectively. Sheath gas and auxiliary gas were set at 30 and 10 arbitrary units, respectively. The capillary temperature was 325 $^{\circ}\text{C}$. The analyzer scanned over a mass range of m/z 81-1 000 for full scan at a mass resolution of 70 000. Data dependent acquisition (DDA) MS/MS experiments were performed with HCD scan. The normalized collision energy was 30 eV. Dynamic exclusion was implemented to remove some unnecessary information in MS/MS spectra.

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1. Chen Y, Zhou Z, Yang W, et al. Development of a data-independent targeted metabolomics method for relative quantification using liquid chromatography coupled with tandem mass spectrometry[J]. Anal Chem, 2017, 89(13): 6954-6962.
2. Dunn W B, Broadhurst D, Begley P, et al. Procedures for large-scale metabolic profiling of serum and plasma using gas chromatography and liquid chromatography coupled to mass spectrometry[J]. Nature Protocols, 2011, 6(7):1060-83.
3. Dunn, Warwick B. et al. Procedures for large-scale metabolic profiling of serum and plasma using gas chromatography and liquid chromatography coupled to mass spectrometry. Nature Protocols. 6, 1060-1083 (2011).
4. Zelena, Eva. et al. Development of a Robust and Repeatable UPLC-MS Method for the Long-Term Metabolomic Study of Human Serum. Analytical Chemistry. 81, 1357-64 (2009).
5. Sangster, T., Major, H., Plumb, R., Wilson, A. J., & Wilson, I. D. A pragmatic and readily implemented quality control strategy for HPLC-MS and GC-MS-based metabonomic analysis. Analyst. 131, 1075-1078 (2006).
6. The HUSERMET Project. Analysis of Serum Samples by UPLC-MS. 2008.
7. Want, E. J. et al. Global metabolic profiling of animal and human tissues via UPLC-MS. Nature Protocols. 8, 17-32 (2013).

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